

TIME MEDICAL SYSTEMS



PICA
the Whole Body MRI System

SYSTEM DATA

PICA

brings a new dimension to MR Imaging.

TIME MEDICAL SYSTEMS

Vision Beyond Imagination.

We See, We Care.



TIME MEDICAL
SYSTEMS

TABLE OF CONTENTS


HTS Technology	5
PICA	6
Specifications	7
Site Planning	9
Room Requirements	11
Typical MRI Floor Plan	12
Chassis & Patient Table	13
Magnetic Fringe Field	14
Performance Highlights	15
Staff Productivity	17
Patient Comfort	19
System Reliability	21
Radio Frequency (RF) Coils	23
MRI Images	25



For the past 30 years, MR technology has relied upon higher magnetic field strengths to achieve a higher Signal-to-Noise Ratio (SNR) in the MR Image. This trend has come to an end, driven mainly by the need for new cost containment strategies and regulatory constraints for imaging above 3 Tesla.

Overcoming the dominant source of noise in an MR Image remains the key challenge. The loss of SNR from the RF coil is large when compared with losses from the sample itself. If RF coils can be optimized to reduce coil loss, they can greatly improve the SNR.

High Temperature Superconducting (HTS) materials and methods have become the ideal approach to improving SNR. The HTS coil is made from YBCO thin film on sapphire substrate. HTS RF Coil technology was developed in the late 1990's by a research team at Columbia University, (New York, USA), led by Professor Q. Y. Ma, PhD.



Phantom image comparison: (a) signal intensity plot through the midline of images of a homogeneous cylindrical phantom. Each curve has been normalized such that the RMS noise is equal in all curves. The vertical lines indicate the range over which the ROI was taken on the images. (b) The corresponding images are shown in the lower panel. Room temperature copper (dotted line), cooled copper (dashed line), and HTS coil (solid line) images are shown from left to right, respectively. The ROI locations for signal and noise are indicated on the middle image.

TIME MEDICAL SYSTEMS was founded by Professor Ma to bring HTS RF Coil technology to the MRI market. This breakthrough innovation led to the era of **low-cost, high-performance MRI Scanners**.


PICA is the first whole body MRI Scanner to take advantage of HTS RF Coil Technology. Operating in conventional mode, or HTS mode, the Hybrid PICA MRI Scanner enables the full range of clinical applications.



The 20 minute Brain MRI Exam is one of the best examples of Hybrid PICA's value, performance & innovation.



HTS Head Coil



Specifications

Magnet System

Type: Permanent, active shims
Material: Neodymium Ferrite Boron (**NdFeB**)
Field Strength: 0.35 Tesla
Opening: 400 mm
Fringe Field:
5 Gauss line 2.5 m axial x 2.5 m radial
See diagrams on page 14.
Size with chassis:
2260 x 2740 x 1996 mm (L x W x H)
Field Stability: < 0.2 ppm/hr
Magnet Weight: 17,200 kg / 37,840 lbs

Gradient System

Amplitude: 25 mT/m (PICA Standard)
Rise Time: 0.50 ms
Slew Rate: 50 T/m/s

** Upgradable to:*

Amplitude: 33 mT/m (PICA Advanced)
Rise Time: 0.37 ms
Slew Rate: 90 T/m/s

RF System

Frequency: 15 MHz \pm 200 KHz
RF Power: 6000 W peak ms, 600 W average
Image BW: 4 KHz – 100 KHz
Preamp Noise Figure < 0.5 dB
Pre-amplifier integrated in all receiving coils
Automatic receiving coil recognition

RF Receive Coils

Quad Head Coil
Quad Neck Coil
Quad Shoulder Coil
Quad Wrist Coil
Quad Body / Spine Coil
Quad Knee Coil
Quad Foot/Ankle Coil

Pulse Sequences

SE 2D & 3D, GRE 2D & 3D, FLASH 2D & 3D, FSE 2D, FSE 3D / HASTE, STIR 2D, IR-FSE 2D (T2 STIR, PD-STIR), FLAIR 2D, Fast-FLAIR 2D, MRA-TOF 2D & 3D, MRV-TOF 2D, 3PT DIXON SE 2D, 3PT DIXON GRE 2D.

Computer System

Host Computer:
Deil™ Optiplex™ with Quad Core processors at 3.1 GHz.
4 GB RAM, 1 TB RAID hard drive
Operating System: Microsoft® Windows® Based
Imaging Software:
Prodiva™ MRI Platform
Data Archive:
DVD-RW (DL) 4.7 GB
Display: 24 inch HD 1080p LCD Monitor
Data Transfer & Handling (HIS/RIS): DICOM 3.0
Printing Support: DICOM Print, All Windows Printer & PDF

Patient Table

Size: 2600 x 800 x 767 mm (L x W x H)
IN/OUT movement – 2000 mm
LEFT/RIGHT movement – 120 mm each side
UP/DOWN movement – 200 mm
Min height from floor – 550 mm
Loading capacity – 200 kg / 440 lbs

Imaging Parameters

FOV: 5 cm – 42 cm, 0.1 mm increments
Slice Thickness:
2D: 1 mm – 10 mm, 0.1 mm increments
3D: 0.25 mm – 10 mm, 0.1 mm increments
Slice Spacing: 0.1 mm increment
Min / Max Phase Matrix: 64 / 512
Min / Max Frequency Matrix: 64 / 512
Bandwidth: 4 kHz – 100 KHz
Rectangular FOV (increment in % in full FOV): 0.1%
Types of Recon Filters:
Ring Filter (On/Off), Histogram Equalization, Image Uniformity

Image Reconstruction

Multithreaded reconstruction
2D > 1000 images/s, 256 x 256
3D > 200 planes/s, 256 x 256 x 64

Image Visualization

Features: Visualization with auto Window & Level, pan, zoom, rotate, flip, multiple image display, ROI, annotation, measurements, cross series marker, color mapping, image analysis, real-time MIP, CINE multi viewport.





Room Requirements





Specification	Equipment Room	Magnet Room	Control Room
Recommended room size (inside dimensions) *	3000 x 3000 mm (9.8 x 9.8 ft)	6500 x 4500 mm (21.3 x 14.7 ft)	3000 x 4000 mm (9.8 x 13.1 ft)
Minimum ceiling height	2800 mm (9.2 ft)	3000 mm (9.8 ft)	no requirement
Floor requirements	antistatic	antistatic, level	antistatic
Total Floor Loading	-	20 kPa 417.7 lbs / sqft	-
Floor Leveling	-	5 mm / 3 m	-
RF shielding	not required	90 dB attenuation over 10 – 100 MHz	not required
Magnetic field shielding	not required	depends on site	not required
Power outlets	3 x 110 Vac or 220 Vac (depends on country) 1 x 380 Vac (3 phase)	3 x 110 Vac or 220 Vac (depends on country)	at least 6 x 110 Vac or 220 Vac (depends on country)
Air conditioning	3 kW	15 kW	3 kW
Humidity	40%-70% without condensation	40%-70% without condensation	40%-70% without condensation
Network outlets	1 (directly connected to Control Room)	No requirement	1 (directly connected to Equipment Room)








* Remarks:


The minimum space requirements vary based on the specific site survey.

Typical MRI Floor Plan











Scan Mode



View Mode



Print Mode





Six Direction Automated Patient Table



UP/DOWN travelling distance 200 mm



IN/OUT travelling distance 2000 mm






LEFT/RIGHT travelling distance \pm 120 mm




Laser Positioning Device \pm 1 mm accuracy

Patient Table	
Cooling	Air
Coolant thermal loading	25 W Ambient (-10°C to +40°C)
Unit weight	450 kg (992 lbs)
Patient loading weight	200 kg (440 lbs)
Shipping weight	600 kg (1322.8 lbs)
Power Input	DC 24V - 40A

System / Power Supply Cabinets	
Cooling	Air cooled with internal fans
Coolant thermal loading	2.8 kW - 5.85 kW
Line voltage	380 Vac 3-phase, neutral $\pm 5\%$, 3N~, 50/60 Hz Optional Power Supply: 400, 420 VAC 3-phase, $\pm 5\%$ 440, 480 VAC 3-phase, $\pm 5\%$
AC power requirements	63 A / 42 kVA
Dimensions (LWH)	910 x 610 x 1735 mm / 670 x 820 x 970 mm
Unit weight	426 kg (939 lbs) / 390 kg (860 lbs)






Item	Specification	Remark
1. Coil type	2 channel, Receive only	Integrated preamp
2. Tuning Method	Passive	
3. Decoupling	Passive & active	From transmission coil
4. Q Factor	>100	F = 15 MHz \pm 200 KHz, loaded
5. Uniformity for Optimal IQ	>90%	
6. Connector	Mixed D type	

Coils	FOV (L x W x H mm)
1. Quad Head Coil	260 x 215 x 250
2. Quad Cervical Spine Coil	180 x 190 x 215
3. Quad Body Coil:	
20" coil	300 x 550 x 325
17" coil	300 x 420 x 280
14" coil	300 x 350 x 240
4. Quad Knee Coil	250 x 180 x 180
5. Quad Foot / Ankle Coil	180 x 115 x 175
6. Quad Wrist Coil	130 x 115 x 70
7. Quad Shoulder Coil	195 x 200 x 215

TIME MEDICAL SYSTEMS




BRAIN T2 FSE Axial 224x256 FOV 24 cm

TIME MEDICAL SYSTEMS



C-SPINE T2 FSE Sagittal 216x256 FOV 26 cm


TIME MEDICAL SYSTEMS



ABDOMEN T1 GRE Axial 160x256 FOV 32 cm

27


TIME MEDICAL SYSTEMS



ABDOMEN T2 FSE Axial 192x256 FOV 32 cm


28

TIME MEDICAL SYSTEMS




T-SPINE T2 IR-FSE Sagittal 160x256 FOV 30 cm

TIME MEDICAL SYSTEMS




L-SPINE T1 SE Sagittal CE 192x256 FOV 26 cm

TIME MEDICAL SYSTEMS



L-SPINE T2 FSE Axial 224x256 FOV 27 cm

TIME MEDICAL SYSTEMS



TIME MEDICAL SYSTEMS

TIME MEDICAL SYSTEMS

USA / Latin America / Europe

Time Medical Systems Inc.
11010 Roselle Street, Suite # 120
San Diego, CA 92121
USA
Tel +1 858 458 0166

Asia Pacific

Time Medical International Pte Ltd.
Block 71, Ayer Rajah Crescent #03-05
Singapore, 139951
Tel + 65 6776 7819

Central / Northern Asia / Africa

Time Medical Limited
G/F, Bio-Informatics Centre
No. 2, Science Park West Avenue
Hong Kong Science Park, Shatin, N.T.
Hong Kong
Tel + 852 2156 1711

China (Shanghai)

Time Medical System (Shanghai) Co., Ltd.
4/F, Block B, 563 Song Tao Road
Zhangjiang Hi-Tech Park
Pudong, Shanghai, 201203
China
Tel + 86 21 5080 4982

China (Beijing)

1226, 12/F, Block A, China Technology Exchange Building
No. 66 North 4th Ring
Haidian, Beijing, 100080
China

China (Southern)

Time Medical Guangdong Co., Ltd.
1st Floor, A2 Building Zone A, #25 Huanshi Avenue South,
Nansha Technology Innovation Center, Nansha, Guangzhou
China, 511458
Tel + 86 20 3905 2216

or contact us at:

sales@time-medical.com
www.time-medical.com

This information is furnished without warranty, representation, or license of any kind, except that it is accurate to the best of Time Medical's knowledge or obtained from sources believed to be accurate. Time Medical does not assume any legal responsibility for use or reliance upon same. Customers are encouraged to conduct their own evaluation. Before using any product, read the product manuals and labels. For products/features under development the information is provided to assist you. Time Medical releases this product for general sale. Time Medical reserves the right to alter prices, specifications, features, capabilities, function, manufacturing release dates, and even general availability of the product at any time.

